

PATIENT & STAFF STATION/CALL POINT

The CarePlus™ multifunction Call Point is a flush wall mount unit with a large push button for various nurse call stations. There are several colours and labels for different levels and alarms, and it comes complete with a reassurance light and tone for patient reassurance.



Technical Data

Socket part number	DAC-W-XXX-Y-ZZ
Ambient temperature	5°C - 55°C
Type of protection	IP 20
Material	PC + ABS
Colour	White
Fixing / mounting	Surface mount on wall or medi duct
Dimension	L 110mm W 70mm D 38mm
Weight	Approx. 80g

Features:

- ▶ Large Easily Identified Nurse Call push button
- ▶ Braille Nurse call button
- ▶ Reassurance light (active call light)
- ▶ Reassurance tone (active call tone)
- ▶ Pull out pendent alarm on socket
- ▶ Stereo 4.5mm jack socket universal input adaptor for normally open connection to any 3rd party device / pendent / bed sensor
- ▶ 2 priorities per call point.
- ▶ Call point button can be one priority (i.e. Call) and jack socket can be alternate priority (i.e. out of bed or staff assist)
- ▶ Raised colourful buttons on polyester Lexan overlay for call and reset
- ▶ No sharp edges
- ▶ Simple installation to surface mount or to wall box
- ▶ Post commissioning surround to reduce over paint and ingress potential
- ▶ Customizable Logo
- ▶ Cut cable alarm

cont. overleaf

AUTHORISED DISTRIBUTOR

USA office:

Europe / UK office:

Australian office:

Variants: Normally Open (momentary closure)

Part number	Description
DAC-M-NC-G	Nurse Call Green + Cancel
DAC-M-NC-R	Nurse Call Red + Cancel
DAC-M-EMG-R	Emergency Red + Cancel
DAC-M-EMG-B	Emergency Blue + Cancel
DAC-M-SA-Y	Staff Assist Yellow + Cancel
DAC-M-NC-G-SA	Nurse Call Green + Cancel with socket
DAC-M-NC-R-SA	Nurse Call Red + Cancel with socket
DAC-M-EMG-R-SA	Emergency Red + Cancel with socket
DAC-M-EMG-B-SA	Emergency Blue + Cancel with socket
DAC-M-SA-Y-SA	Staff Assist Yellow + Cancel with socket

Pin No.	Name	Description
1	Power	+24v
2	Reassurance	Reassurance LED on call point
3	NC	
4	NC	
5	NC	
6	NC	
7	Input	Patient call
8	Input	Call Cancel

The CarePlus™ Call Point is UL1069 and DIN 0834 certified, and meets AS3811.



AUTHORISED DISTRIBUTOR

USA office:

1000 N West Street, Suite 1200
Wilmington
Delaware 19801
t: +1 302-295-4845
e: info.usa@niq-health.com

Europe / UK office:

16 Upper Woburn Place
London WC1H 0AF
t: +44 (0)20 7193 0072
e: info.euorpe@niq-health.com

Australian office:

Unit 2 / 3 King Edward Road
Osborne Park, Perth
Western Australia 6017
t: +61 8 9242 7407
e: info.aus@niq-health.com